TGR-Y MR 001-4

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HARMONY OF LIFE

DEAR CUSTOMER,

Thank you for purchasing our product. We hope that the equipment you have selected will meet your expectations. We would like to share some important information with you.

Timago International Group is a Polish company which distributes medical, rehabilitation and orthopaedic products to pharmacies and medical stores.

Our mission is to ensure your comfort by offering safe, functional, high-quality products.

We have many years of experience in this sector. Timago products reach a wide range of users in Poland and throughout Europe through our distributors (medical retailers).

Timago's company culture promotes openness to your needs. Our goal is to promote safe and correct use of our products for treatment and rehabilitation in order to improve comfort and quality of life.

We invite you to access our website for more information: **timago.com**.

1. INTRODUCTION

Prior to use of our products you should inspect the equipment for any defects which may have occurred during transportation. If this happens, please contact your retailer.

Next, please read these instructions which include important user information.

Please contact us if you have any questions concerning the use of the equipment. Our advisors will provide you with any information you require on an ongoing basis.

Please note that if you have questions other than those concerning the use of our products, you should contact your doctor, nurse or physical therapist.

2. PRODUCT DESCRIPTION

The anti-decubitus strip mattress with pump is a set which includes a pressure alternating mattress and a compressor (air pump). The mattress is made of durable PVC in the form of transverse chambers which are alternately inflated in cycles by the compressor (air pump). By changing support points of the person lying on the bed, this set minimizes pressure on the parts of the body which are the most at risk from bedsores and makes it possible to avoid or reduce skin damage. The PVC used in the mattress is made of a non-toxic and hypoallergenic material.

The set includes:

- Soft Wind VI alternating pressure pump,
- Tube chamber mattress,
- Air hoses,

- Spare chamber,

3. INDICATIONS

Anti-decubitus strip mattresses with pump are aids for people at risk of bedsores. They are used for people who are bedridden for an extended period of time, in whom peripheral circulatory system disorders occur, and also for individuals with signs of subcutaneous or open skin damage caused by long-term pressure. These mattresses are used to prevent bedsores and treat stage I (according to NPUAP / EPUAP) bedsores. The mattress works by cyclically changing support points of the person lying on the bed, minimizing pressure on the parts of the body which have the greatest risk of bedsores. This improves circulation of the tissues of the body parts which remain in contact with the bed surface for the longest period. It should be noted that the mattress operating principle described above cannot eliminate all causes of skin lesions; however, it can significantly limit their development.

Contraindications for use of alternating pressure mattresses are determined by a physician based on a physical examination and medical history. As a general guideline, the mattress is contraindicated if alternating changes of support points cause pain and an exacerbation of the patient's condition.

4. ASSEMBLY AND PRIMARY INDICATIONS FOR USE

FOR THE ASSEMBLY AND SET-UP DESCRIBED BELOW, THE USER MAY REQUIRE THE HELP OF ANOTHER PERSON. YOU SHOULD ASK FOR HELP FROM A PERSON WITH PROPER SKILLS WHO SHOULD FOLLOW INFORMATION INCLUDED IN THE INSTRUCTIONS.

THE MATTRESS SHOULD NOT BE USED IF IT IS NOT FULLY INFLATED.

STEP 1

Place the mattress on the bed directly on top of the regular mattress. The flaps at each end of the pressure relief mattress should be tucked under the regular mattress to limit movement of one mattress relative to the other. It is recommended that the mattress be secured with a hygienic vapour-permeable mattress pad or cotton sheet in order to prevent it from being soiled.

STEP 2

The hook of the pump should be opened to the desired width and hung on the bed frame.

STEP 3

The pump and the mattress should be connected with air hoses so that they are not kinked or trapped under the mattress. Make sure that the hoses are properly fastened and do not come out from the ports.

STEP 4

Plug in the pump using the power cord, taking extra care while doing this.

STEP 5

To start the pump, press the ON button on the control panel. The pressure relief mattress will then start to inflate.

STEP 6

Set the pressure dial at the maximum value. The preparation time (the first inflation) prior to use takes 24 minutes (two full cycles).

YOU MUST WAIT 24 MINUTES FROM STARTING THE MATTRESS BEFORE USE.

STEP 7

Following complete inflation, the pressure dial should be set at a value appropriate to the patient's weight. The pressure relief mattress with air pump will then be ready to use.

THE PRESSURE SHOULD BE APPROPRIATE FOR THE USER'S INDIVIDUAL NEEDS DETERMINED BY HIS/HER GENERAL CONDITION AND COMPLAINTS ASSOCIATED WITH THEIR CONDITION. TO OBTAIN THIS INFORMATION YOU SHOULD CONTACT YOUR DOCTOR.

THE AIR PUMP OPERATES BY ALTERNATELY PUMPING AIR TO MATTRESS AIR CELLS AND FOR THIS REASON NEEDS TO BE ON ALL THE TIME. THE OFF BUTTON ON THE CONTROL PANEL IS ONLY FOR TURNING OFF THE PUMP WHEN THE MATTRESS IS NOT IN USE. BEFORE USING THE MATTRESS AGAIN, THE STARTING PROCEDURE MUST BE PERFORMED (SEE STEPS 6 AND 7).

5. TECHNICAL SPECIFICATIONS

	MATTRESS	AIR PUMP
MATTRESS DIMENSIONS	200 X 86 X 9,5 CM	
WEIGHT	3,6 KG	1,2 KG
MAXIMUM WEIGHT CAPACITY	130 KG	
CYCLE DURATION	12 MIN	
MODEL	MADE OF PVC IN THE FORM OF TRANSVERSE CHAMBERS	SOFT WIND VI
PRESSURE RANGE		40-110 MMHG
FLOW		6-7 L/MIN
POWER		AC 220V 50/60 HZ 0,1A
FUSE		T1A 250V
WATTAGE		8W
OPERATING ENVIRONMENT	TEMPERATURE: +10°C-40°C, HUMIDITY: 10%-90%	
TYPE OF OPERATION	CONTINUOUS - CYCLICAL CONTINUOUS, ALTERNATING AIR CELL INFLATION	

ELECTRICAL DIAGRAM



DIAGRAM OF THE SOFT WIND VI AIR PUMP

6. SAFETY RULES

BEFORE USING THE MATTRESS YOU MUST TALK TO YOUR DOCTOR IN ORDER TO DETERMINE ANY CONTRAINDICATIONS FOR USE. THE MATTRESS SHOULD NOT BE USED WITHOUT CONSULTING A DOCTOR. THIS IS DUE TO THE RISK OF ADVERSE EFFECTS, IN PARTICULAR:

- IN PATIENTS ON CERVICAL TRACTION**
- WHERE PATIENTS HAVE MULTILEVEL INSTABILITY OF A GIVEN SEGMENT OF THE SPINE**

IN PATIENTS WITH A PACEMAKER, THE PUMP SHOULD BE PLACED AS FAR FROM THE CHEST AS POSSIBLE.

ONE MATTRESS SHOULD ONLY BE USED FOR ONE PATIENT.

DO NOT EXCEED THE MAXIMUM WEIGHT CAPACITY (130 KG).

NO REHABILITATION PROCEDURES SHOULD BE PERFORMED ON THE MATTRESS.

OTHER PEOPLE SHOULD NOT SIT ON THE MATTRESS.

THE MATTRESS MAY ONLY BE USED ON FLAT SURFACES. THE MATTRESS SHOULD BE PROTECTED AGAINST SHARP OBJECTS AND SHOULD NOT BE LAID DIRECTLY ON BED SPRINGS OR BED FRAMES WITH SHARP EDGES.

THE AIR PUMP SHOULD NOT BE USED IF THE POWER CORD OR PLUG IS DAMAGED, OR IN CASE OF APPARENT PRODUCT FAILURE, FOR EXAMPLE, INSUFFICIENT PUMPING FORCE. THE AIR PUMP SHOULD NOT BE CONNECTED TO DAMAGED EXTENSION CORDS OR ELECTRIC OUTLETS.

KEEP THE MATTRESS FROM TOUCHING HEAT SOURCES AND DO NOT PLACE THE AIR PUMP CLOSE TO HEAT SOURCES. KEEP AWAY FROM OPEN FLAMES AND FLAMMABLE SUBSTANCES.

THE PUMP SHOULD NOT BE USED IN AN IMPROPER OPERATING ENVIRONMENT OR COME IN CONTACT WITH WATER OR OTHER FLUIDS. THE AIR PUMP SHOULD NOT BE PLACED UNDER THE BED OR IN PLACES WHERE THE AIR FLOW TO THE AIR PUMP MAY BE LIMITED, E.G. UNDER BEDDING. THE AIR PUMP COVER MUST NOT BE REMOVED DUE TO THE RISK OF ELECTRIC SHOCK.

IF THE MATTRESS IS NOT IN USE, IT IS IMPERATIVE TO TURN OFF THE AIR PUMP USING THE OFF BUTTON ON THE CONTROL PANEL; AFTER THE AIR PUMP IS TURNED BACK ON, THE OPERATING INSTRUCTIONS MUST BE FOLLOWED.

DO NOT PULL ON THE HOSES CONNECTING THE AIR PUMP WITH THE MATTRESS. THIS CAN CAUSE DAMAGE TO THE MATTRESS.

DO NOT ATTEMPT TO REPAIR OR RETROFIT THE MATTRESS ON YOUR OWN.

WE ONLY RECOMMEND USING REPLACEMENT PARTS FROM THE MANUFACTURER. PLEASE CONTACT YOUR RETAILER AS REQUIRED.

A WORN OUT OR DAMAGED AIR PUMP SHOULD BE TAKEN TO AN ELECTRICAL OR ELECTRONIC RECYCLING CENTRE. DO NOT DISPOSE OF THE AIR PUMP WITH OTHER MUNICIPAL WASTE.

TIMAGO INTERNATIONAL GROUP IS NOT RESPONSIBLE FOR THE CONSEQUENCES OF IMPROPER PRODUCT MAINTENANCE, FAILURE TO FOLLOW SAFETY RULES OR USE INCONSISTENT WITH ITS PURPOSE.

7. MAINTENANCE

You should personally check the condition of the air pump and the mattress regularly.

- Prior to each use: check whether air hoses are properly fastened in the ports, that air is not leaking from the mattress and that the air pump is working smoothly and efficiently.
- Every 6 months: It is recommended to have an inspection performed by a service facility. You should contact the manufacturer or your retailer as required.

CLEANING

The mattress and the air pump should be cleaned with a damp cloth using mild cleaning products certified by the National Institute of Hygiene, paying attention not to get it wet, and wiped dry after each cleaning. The use of strong cleaning products or detergents can damage the product.

NO HEAT SOURCES SHOULD BE USED FOR DRYING.

PRIOR TO CLEANING THE AIR PUMP, YOU MUST TURN OFF THE AIR PUMP USING THE OFF BUTTON ON THE CONTROL PANEL AND UNPLUG THE POWER CORD.

8. WARRANTY INFORMATION

All products distributed by our company are under a warranty, the terms of which are described in the warranty card received at the time of purchase. Please keep in mind that for warranty purposes, proof of purchase (receipt or invoice) should be retained.

As a user friendly company, we only supply products which have been thoroughly tested with respect to materials, quality and functionality.

Please contact us if you have any questions concerning our servicing procedure.

9. LABELS

Explanation of terms found in product labels:

10. RECOMMENDED PRODUCTS

Products offered by Timago International Group include:

- Orthotics and joint stabilizers
- Manual and special-use wheel-chairs
- Walkers and walking aids
- Crutches and canes
- Rehabilitation beds
- Adaptive and bathroom equipment
- Equipment and accessories to prevent bedsores
- Stabilizing and immobilizing belts
- Electro-medical equipment
- Disposable medical equipment

Our complete line of products is available from our distributors.

Learn more at: timago.com.

	REFERENCE NUMBER		MANUFACTURER
	LOT NUMBER		MANUFACTURE DATE
	SERIAL NUMBER		MEDICAL DEVICE
	NOTE		PLEASE READ THE INSTRUCTIONS
	THE MANUFACTURER CHECKED FOR COMPLIANCE WITH ESSENTIAL MEDICAL DEVICES REQUIREMENTS.		

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