

- EN** Folding wheelchair in painted steel
F Fauteuil roulant pliant en acier vernis
ES Silla de ruedas plegable en acero pintado
P Cadeira de rodas dobrável em aço pintado

Kiddy

Single cross frame - Children

Cod. **CP880-35**



Height adjustable handle
86÷99 cm



Supplied with



Elevating footrests

Optional accessories



Cod. **CPA600**
Adjustable IV Pole



Cod. **CPA605**
Chromed steel oxygen
bottle holder Ø Max
oxygen bottle 12 cm



Cod. **CPA606**
Polyester Oxygen bottle
holder - Ø oxygen bottle
Max 12 cm



Cod. **CPA507**
Table for wheelchair
58x38 cm



Cod. **CPA501**
Linear table
55x32 cm



Cod. **RP247**
Safety belt with
spring buckle

Spare parts



Cod. **CPR810**
Flip up armrest -PAIR



Cod. **CPR821**
PU rear wheel Ø 56
cm - SINGLE



Cod. **CPR820**
PVC front wheel Ø
15 cm - SINGLE



Cod. **CPR827**
Armrest padding only



Cod. **CPR815**
Lateral footrests



Cod. **CPR817**
Elevating footrests

Technical features

SPECIFICATIONS

Front wheels	PVC, black, painted steel fork														
Rear wheels	Plastic radius, grey PU wheels														
Upholstery	Padded nylon- Fire reaction tested, according to EN1021-1/2														

WEIGHT AND DIMENSION

Measures in cm															
CP880-35	94	74	36	86÷99	62	42	37÷42	35	37	56	31	14,5	56	15	- 100