

## INSTRUCTIONS/ ΟΔΗΓΙΕΣ

Ref.: 0808911

### FOAM MATTRESS 3-SECTION)/

### ΣΤΡΩΜΑ ΠΟΛΥΣΠΑΣΤΟ ECONOMIC ΑΦΡΟΛΕΞ

#### EN

The mattress has been compressed to reduce its volume during transportation.

It needs 24 hours out of packing to come back to its original dimensions.

(The Mattress is not Washable)

The Washing Instructions Refer to the Cover ONLY:

- Hands Wash Only, using mild soap.
- Do Not Bleach
- Do Not Iron.
- Do Not Tumble Dry.
- Do Not Dry Clean.
- Do Not Dry Wring.

#### GR

Το στρώμα έχει συμπιεστεί για τη μείωση του όγκου του κατά τη μεταφορά.

Χρειάζεται 24ώρες εκτός συσκευασίας για να επανέλθει στις αρχικές του διαστάσεις.

Οι Οδηγίες Πλυσίματος Αφορούν μόνο το Κάλυμμα:

- Πλύσιμο μόνο στο χέρι, χρησιμοποιώντας ήπιο απορρυπαντικό.
- Απαγορεύεται το χλώριο.
- Απαγορεύεται το σιδέρωμα.
- Μην Χρησιμοποιείτε στεγνωτήριο.
- Απαγορεύεται το στεγνό καθάρισμα.
- Απαγορεύεται το στύψιμο.



**WORLD HEADQUARTERS:**  
AKROTIRI CHANIA CRETE P.C 73100  
T: +30 28210 63222

**DISTRIBUTION CENTER: CENTRAL, SOUTHERN &  
INSLAND GREECE:**  
SXINEZA ASPROPIRGOS, Attica P.C 19300  
T: +30 210 4614006

**DISTRIBUTION CENTER: NORTHERN GREECE & BALKAN  
AREA:**  
8 ARKADIOU STR. KALOCHORI, THESSALONIKI P.C 57009  
T: +30 2310 785668

**DISTRIBUTION CENTER: CENTRAL & EASTERN CRETE:**  
P' STR. INDUSTRIAL ZONE OF HERAKLION  
CRETE, P.C 71601, T: +30 2810 333151

